

# PRATIKO AUTO BRAKE SELECTION FORM

To ease the ordering process, please use one page per braking system you order  
**Select the configuration of the wheelchair to equip (put an "X" in the appropriate case)**

## 1. Manufacturer and Model (KNOWN COMPATIBILITY)

X	Manufacturer	Model	installation kit #	N	X	Manufacturer	Model	installation kit #	N
	KI MOBILITY	CATALYST 4C	PRA-U-CAT4C			MOTION COMPOSITE	CHRONOS 1	PRA-U-CHKIT*COBALT	
	KI MOBILITY	CATALYST 5TI	PRA-U-CHKIT*MYON			MOTION COMPOSITE	CHRONOS 2	PRA-U-CHKIT*COBALT	
	KI MOBILITY	CATALYST 5Vx	PRA-U-CHKIT*PRIMA	2		MOTION COMPOSITE	COBALT 1	PRA-U-CHKIT*COBALT	
	KI MOBILITY	CATALYST 5C	PRA-U-CHKIT*MYON			MOTION COMPOSITE	COBALT 2	PRA-U-CHKIT*COBALT	
	EVEREST COPY	FIXED AXLES	PRA-U-CHKIT*FXAXL			MOTION COMPOSITE	HÉLIO XX (ALL CONFIG))	PRA-U-CHKIT*HELIOS	
	DRIVE	FIXED AXLES	PRA-U-CHKIT*FXAXL			MOTION COMPOSITE	MOVE	PRA-U-CHKIT*COBALT	
	EVEREST ET JENNINGS	FIX AXLES (L3, L4)	PRA-U-CHKIT*FXAXL			MOTION COMPOSITE	PLATINE 1	PRA-U-CHKIT*HELIOS	
	EVEREST ET JENNINGS	METRO	PRA-U-CHKIT*PRIMA	2		MOTION COMPOSITE	PLATINE 2	PRA-U-CHKIT*HELIOS	
	EVEREST ET JENNINGS	P2+	PRA-U-CHKIT*P2+	1		ORTHO FAB	CH	PRA-U-CHKIT*CH	
	GORDEN	ALL FIXED AXLES	PRA-U-CHKIT*FXAXL			ORTHO FAB	PRIMA	PRA-U-CHKIT*PRIMA	2
	GUARDIAN	EASYCARE 2000	PRA-U-CHKIT*P2+	1		ORTHO FAB	PRIMA 4	PRA-U-CHKIT*PRIMA	2
	GUARDIAN	EASYCARE 4000	PRA-U-CHKIT*P2+	1		PHYSIPRO	XL5	PRA-U-CHKIT*PRIMA	2
	INVACARE	9000 XT, SL, L	PRA-U-CHKIT*FXAXL			PRIDE MOBILITY	LITE STREAM	PRA-U-CHKIT*PRIMA	2
	INVACARE	ACTION	PRA-U-CHKIT*ACTIONINV			PRIDE MOBILITY	QUANTUM REHAB	PRA-U-CHKIT*PRIMA	2
	INVACARE	MVP	PRA-U-CHKIT*MYON			PRIDE MOBILITY	STYLLUS	PRA-U-CHKIT*FXAXL	
	INVACARE	MYON	PRA-U-CHKIT*PRIMA	2		SUNRISE/BREEZY	BREEZE 600	PRA-U-CHKIT*BREEZE	
	INVACARE	PATRIOT	PRA-U-CHKIT*PATRIOT			SUNRISE/BREEZY	QUICKIE I ET II	PRA-U-CHKIT*PRIMA	2
	INVACARE	X4	PRA-U-CHKIT*PRIMA	2		SUNRISE/BREEZY	QS	PRA-U-CHKIT*PRIMA	2
	INVACARE	TRACER SX5	PRA-U-CHKIT*FXAXL			SUNRISE/BREEZY	ULTRA 4	PRA-U-CHKIT*ULT4	
	MEYRA	EUROCHAIR	PRA-U-CHKIT*EURO			SUNRISE/BREEZY	QX	PRA-U-CHKIT*MYON	
	MAPLE LEAF	ACCESS	PRA-U-CHKIT*FXAXL			SUNRISE/BREEZY	LXI/LX	PRA-U-CHKIT*PRIMA	
	MAPLE LEAF	EZRIDE	PRA-U-CHKIT*PRIMA	2					
	MAPLE LEAF	NRG+	PRA-U-CHKIT*FXAXL						
	MAPLE LEAF	NRG+ RECLINER	PRA-U-CHKIT*PRIMAMOD	3					
	MAPLE LEAF	NRG+ GOLD	PRA-U-CHKIT*PRIMAMOD	3					
	MAPLE LEAF	SWIFT	PRA-U-CHKIT*PRIMA	2					

- NOTES:  
 1, P2+ KIT CAN BE REPLACE WITH FXAXL KIT  
 2, PRIMA KIT CAN BE REPLACE WITH MYON KIT (BUT THE MYON KIT CANNOT BE REPLACE WITH PRIMA KIT)  
 3, PRIMA KIT MODIFIED TO ATTACH AT THE HEMI (LOWEST ADJUSTMENT) POSITION

If you cannot find the manufacturer or the model of the wheelchair you wish to install please contact us at (819) 826-3300. Most of the wheelchair fall under one of listed kit we will be please to indicate the one to choose and to update our list

## 1B. LIST OF INCOMPATIBLE WHEELCHAIRS Manufacturer and Model

Manufacturer	Model	DETAILS
SUNRISE/BREEZE	QXI (TYPE 3)	INCOMPATIBLE AT THE AXLE ATTACHMENT

## 2. Wheels size (usually the wheelsize is indicated on the side of the tire). Otherwise measure the hight of the tire from the ground to the top.

X	Wheel Diameter (in)	blocking disk kit #
	18	PRA-U-CH018
	20	PRA-U-CH020
	22	PRA-U-CH022
	24	PRA-U-CH024

## 3. Hand rims attachment to the wheel (type)

Hollow ties



Tab type



Need the installation bracket: PRA-U-KITLANG

## 4. Identification

Your Name: \_\_\_\_\_

Reference (Mr. Smith room....): \_\_\_\_\_

how can we reach you

Phone: \_\_\_\_\_

email: \_\_\_\_\_

shipping adress:

Name \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_

Province /state: \_\_\_\_\_

Country \_\_\_\_\_

Postal Code (ZIP) \_\_\_\_\_