

PRATIKO AUTO BRAKE SELECTION FORM

To ease the ordering process, please use one page per braking system you order
 Select the configuration of the wheelchair to equip (put an "X" in the appropriate case)

Take note that the HCPCS CODES provided are our interpretation of the codes difinitions. The use of HCPCS does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy

1. COMMON PART INCLUDED IN A KIT

DETECT SYST ASSBLY	PRA-U-CHBASE	E2228
--------------------	--------------	-------

2. Manufacturer and Model (KNOWN COMPATIBILITY)

X	Manufacturer	Model	Installation kit #	HCPC
	KI MOBILITY	CATALYST 4C	PRA-U-CAT4C	E2228
	KI MOBILITY	CATALYST 5Ti	PRA-U-CHKIT*PRIMA	E2228
	KI MOBILITY	CATALYST 5Vx	PRA-U-CHKIT*PRIMA	E2228
	KI MOBILITY	CATALYST 5C	PRA-U-CHKIT*PRIMA	E2228
	EVEREST COPY	FIXED AXLES	PRA-U-CHKIT*FXAXL	E2228
	DRIVE	FIXED AXLES	PRA-U-CHKIT*FXAXL	E2228
	EVEREST ET JENNINGS	FIX AXLES (L3, L4)	PRA-U-CHKIT*FXAXL	E2228
	EVEREST ET JENNINGS	METRO	PRA-U-CHKIT*PRIMA	E2228
	EVEREST ET JENNINGS	P2+	PRA-U-CHKIT*FXAXL	E2228
	GORDEN	ALL FIXED AXLES	PRA-U-CHKIT*FXAXL	E2228
	GUARDIAN	EASycARE 2000	PRA-U-CHKIT*FXAXL	E2228
	GUARDIAN	EASycARE 4000	PRA-U-CHKIT*FXAXL	E2228
	INVACARE	9000 XT, SL, L	PRA-U-CHKIT*FXAXL	E2228
	INVACARE	ACTION	PRA-U-CHKIT*PRIMA	E2228
	INVACARE	MVP	PRA-U-CHKIT*PRIMA	E2228
	INVACARE	MYON	PRA-U-CHKIT*PRIMA	E2228
	INVACARE	PATRIOT	PRA-U-CHKIT*PATRIOT	E2228
	INVACARE	X4	PRA-U-CHKIT*PRIMA	E2228
	INVACARE	TRACER SX5	PRA-U-CHKIT*FXAXL	E2228
	MEYRA	EUROCHAIR	PRA-U-CHKIT*EURO	E2228
	MAPLE LEAF	ACCESS	PRA-U-CHKIT*FXAXL	E2228
	MAPLE LEAF	ERIDE	PRA-U-CHKIT*PRIMA	E2228
	MAPLE LEAF	NRG+	PRA-U-CHKIT*FXAXL	E2228
	MAPLE LEAF	NRG+ RECLINER	PRA-U-CHKIT*PRIMA	E2228
	MAPLE LEAF	NRG+ GOLD	PRA-U-CHKIT*PRIMA	E2228
	MAPLE LEAF	SWIFT	PRA-U-CHKIT*PRIMA	E2228

NOTES:

If you cannot find the manufacturer or the model of the wheelchair you wish to install please contact us at (819) 839-3303. Most of the wheelchair fall under one of listed kit we will be please to indicate the one to choose and to update our list



2B. LIST OF INCOMPATIBLE WHEELCHAIRS Manufacturer and Model

Manufacturer	Model	DETAILS

3. Wheels size (usually the wheelsize is indicated on the side of the tire). Otherwise measure the hight of the tire from the ground to the top.

X	Wheel Diameter (in)	HCPC	blocking discs kit #
	18	E2228	PRA-U-CH018
	20	E2228	PRA-U-CH020
	22	E2228	PRA-U-CH022
	24	E2228	PRA-U-CH024
	26	E2228	PRA-U-CH026

4. Hand rims attachment to the wheel (type)



Need the installation bracket: PRA-U-KITLANG HCPC E2228

5. Identification

Your Name: _____

Reference (NAME TAG / WORK ORDER #...): _____

how can we reach you

Phone: _____

email: _____

shipping address:

Name: _____

Address: _____

City: _____

Province /state: _____

Country: _____

Postal Code (ZIP): _____